Pinbohrgerät M17 compact

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Pin drilling unit M17 compact

 Instruction manual

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Warning: Check that all components are complete and undamaged before use!

Technical specifications

- Dimensions D/W/H: 190 x 380 x 270 mm
- Weight: 6.5 kg
- Connected loads: 230 V/0.32 W
- Rotary speed: 3,600 RPM
- Laser class: 3 A < 5 mW

Package contents

- **DIVARIO®** pin drilling unit M17 compact (art. no. 14500-C)
- **DIVARIO®** carbide burr, 1 pc (art. no. 10171)
- **DIVARIO®** Duo plate support (art. no. 14083)
- Impression carrier (art. no. 14090)
- Assembly tools

Assembly tools

- Never change the burr when the machine is switched on
- Keep fingers and hands away from the area of the laser dot when the machine is switched on – danger of injury!
- The pin drilling unit M17 compact is not designed for use as follows:
 - In explosion hazard areas
 - For medical applications
 - For working with damp or wet materials
- The Mälzer Dental product warranty will be voided by any handling of the pin drilling unit M17 compact other than that specified in the user manual
- Observe the regulations of the relevant employer's liability insurance association when using the M17 compact pin drilling unit:
 - Always wear protective goggles
 - Ensure adequate workplace lighting
 - Use an extraction unit
- **Warning:** All repairs or other work on the machine must be carried out solely by Mälzer Dental or qualified personnel authorised by Mälzer Dental!

General information & care instructions

- Before commissioning, check that the mains power specifications match those on the rating plate.
- Clean the pin drilling unit from time to time with a soft cloth. Do not use compressed air or detergents.
- If the machine develops a fault, always return the unit directly to Mälzer Dental or an authorised dealer.
- High-quality technical machinery of this kind requires an appropriate repair service and specially trained personnel. Mälzer Dental guarantees that repairs will be performed properly with original replacement parts and offers a six-month guarantee on all repairs!

Legal information

The information presented here is correct to the best of our knowledge. Errors and omissions excepted. We accept no legal responsibility for the consequences of inaccurate statements or improper usage.

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Instruction manual

Commissioning

Attention – do not lift the pin drilling unit by grasping the laser arm but solely under the magnet plate!

Electrical connection, incl. power plug (220 V/320 W) and main switch used to switch the pin drilling unit off and on.

Release button for starting the drilling procedure

Setting the drilling depth

Note: the optimum drilling depth is factory-set!
The drilling depth is regulated by using an Allen key to rotate the set screw located under the burr (centre of the axis) to the chosen position.

Clockwise rotation: Screw is lowered = drilled hole is narrower
Anticlockwise rotation: Screw is raised = drilled hole is wider

Changing the burr

Switch off and unplug the machine. Using the Allen key provided, disengage the two grub screws to the side of the brass spindle and pull the burr upwards to remove.

Now insert the new burr from above and retighten the two grub screws. To set the drilling depth, proceed as shown in Figure 5.

Switch the machine back on and carry out a test drill.


Placing and blocking-out the impression + drill procedure

Trim the impression, blocking-out missing areas with wax or silicone.


Lie the positioning plate parallel to the impression.

Use putty/modelling material to fix the impression on the impression carrier.

To orient the impression, use the small or large positioning plate.

Block-out missing areas with putty, so as to avoid spillage of plaster mush when pouring out.


Clamp the base plate into the pin drilling unit's plate support.

Clamped base plate.

Place the impression carrier on the plate support.

Place the impression carrier incl. plate support on the pin drilling unit's magnetic plate.


Set/target the desired drilling position using the laser pointer.

WARNING: The laser beam must be positioned within the **DIVARIO®** base plate!

Press the button to start drilling – the plate support is held in place electromagnetically.

WARNING: If the drill position/laser beam is outside the dental arch, then the burr may break, due to coming into contact with the plate support or the retention disc!

Base plate with drilling complete.

Populate the base plate with pins.

Once pinned, check the base plate by placing it in the impression.

Pouring out and finishing

Fill the impression with plaster mush. To improve flow properties, increase shaker vibration strength as required.


To avoid trapping air bubbles when dipping the base plate, fill the impression with more plaster than necessary.

Invert the base plate onto the guide in the carrier plate and lower gently.


Let plaster set slightly before carefully remove excess around base plate.

After 25 minutes, carefully remove the impression incl. base plate.


Demoulding, trimming and sawing

Tap carefully to free the dental arch from the base plate.


As an alternative to the plaster trimmer, work the dental arch freehand (by using the finishing strips, order no.: 15100-120).

Mark sawn sections beforehand for better orientation.

Saw the dental arch into sections with a diamond blade.

Work the individual segments with a finishing strip or plaster bur and clean to remove plaster dust.


Completed **DIVARIO**® sawn model, seated perfectly on the base plate.

Position the **DIVARIO**® secondary plate before mounting, so as to guarantee a perfect fit.

Finished, mounted work in the articulator.

For large segments, mark additional relief cuts (see sketch). Pins are represented by red dots in the sketch.

Schlesierweg 27 · D-31515 Wunstorf (Germany)
Phone: +49 (0) 50 33 / 963099-0 · Fax: +49 (0) 50 33 / 963099-99
E-Mail: info@maelzer-dental.de · Internet: www.maelzer-dental.de

maelzer dental 